



March 25, 2024

Travel Compliance Letter for ResMed Therapy Devices

To whom it may concern, ResMed respiratory assistive devices (ventilators, respirators, and CPAP machines) are medical devices prescribed for use on or by the order of a physician and are cleared for sale by the U.S. Food and Drug Administration, a branch of the Federal Government of the United States of America, and the Australian Therapeutics Goods Administration.

The devices do not contain any parts or materials known to cause or interfere with the safe operation of a commercial aircraft.

The U.S. Department of Transportation (DOT) Final Rule, “Non-discrimination on the Basis of Disability in Air Travel” (73 FR 27614 which updates Title 14 CFR Part 382), effective May 13, 2009, provides important new requirements for the accommodation of passengers with respiratory assistive devices.

Specifically, respiratory assistive devices may be used on board an aircraft, without further testing by the carrier, provided they have been tested for electromagnetic compatibility (EMC) in accordance with the current version of RTCA/DO-160, Section 21, Category M.

ResMed has successfully completed testing for the respiratory assistive devices described in this letter.


The devices listed below comply with RTCA/DO-160G, Section 21, Category M, and are considered FAA compliant. *This information may be engraved on the device or a sticker label applied to the device; either method for identifying FAA compliance is acceptable and in accordance with DOT regulations.*

ResMed product	Product codes
AirCurve™ 10 ASV	37042, 37043, 37044, 37092, 37215, 37216, 37229, 37230, 37291, 37410, 37411, 37384, 37144
AirCurve 10 CS PaceWave™	37060, 37061, 37062, 37069, 37070, 37114, 37275, 37276, 37277, 37278, 37279, 37280, 37281, 37282, 37283, 37284, 37341, 37354, 37369, 37438, 37459, 28517, 28516
AirCurve 11 ASV AirCurve 11	39012, 39013, 39023
VAuto AirCurve 11 ASV	39018, 39019, 39020, 39144, 39164, 39182, 39183, 39500, 39502, 39505
PaceWave AirCurve 11 S	39146, 39147, 39168, 39184, 39185, 39193, 39501, 39507, 39506, 39508
AirCurve 11 ST Aircurve 10	39165, 39503
CS-A Aircurve 10 CS-A	39166, 39504
PaceWave™	28281
Aircurve 10 CS-A TJ	28282, 28283, 28284
AirCurve 10 S	28280
	37045, 37046, 37047, 37048, 37068, 37093, 37213, 37214, 37227, 37228, 37250, 37251, 37290, 37340, 37351, 37397, 37408, 37409, 37116, 37107
AirCurve 10 ST	37056, 37057, 37058, 37073, 37095, 37306, 37307, 37348, 37349, 37373, 37374, 37375, 37412, 37413, 37385, 37145, 37111
AirCurve 10 ST-A	28210, 28211, 28219, 28323, 28326, 28285, 28340

ResMed product	Product codes
AirCurve 10 VAuto	37050, 37051, 37052, 37053, 37054, 37055, 37059, 37094, 37109, 37211, 37212, 37225, 37226, 37238, 37239, 37248, 37249, 37261, 37262, 37267, 37272, 37289, 37391, 37392, 37406, 37407, 37450, 37481, 37383, 37143, 37110, 37108, 37164
AirMini™ AutoSet	38112, 38113, 38119, 38121, 38122, 38139, 38140, 38145, 38151
AirSense™ 10 AutoSet™	37027, 37028, 37029, 37030, 37031, 37032, 37033, 37064, 37065, 37067, 37090, 37101, 37102, 37207, 37208, 37221, 37222, 37236, 37237, 37245, 37246, 37254, 37258, 37259, 37266, 37271, 37273, 37274, 37287, 37294, 37315, 37316, 37319, 37321, 37327, 37339, 37352, 37367, 37372, 37402, 37403, 37415, 37436, 37445, 37448, 37453, 37457, 37382, 37142, 37497, 37496, 37122, 37503, 37146, 37147, 37504, 37505, 37147, 37103, 37106, 37126, 37160, 37493, 37426, 37127, 37128
AirSense 11 AutoSet	39000, 39001, 39007, 39140, 39141, 39163, 39180, 39181, 39190, 39191, 39420, 39463, 39466, 39480, 39483, 39485, 39488, 39517, 39105, 39108,
AirSense 11 Auto	39520
AirSense 10 AutoSet for Her	39315 37034, 37035, 37036, 37037, 37039, 37040, 37041, 37066, 37091, 37104, 37105, 37209, 37210, 37223, 37224, 37247, 37260, 37288, 37295, 37317, 37318, 37320, 37322, 37328, 37353, 37368, 37390, 37396, 37404, 37405, 37420, 37421, 37422, 37423, 37430, 37434, 37437, 37446, 37449, 37454, 37458, 37494, 37129, 37130, 37131
AirSense 10 CPAP	37014, 37015, 37016, 37017, 37018, 37063, 37072, 37088, 37203, 37204, 37217, 37218, 37285, 37311, 37312, 37325, 37395, 37398, 37399, 37440, 37461, 37497
AirSense 11 CPAP	37461, 37497
AirSense 10 Elite	39002, 39003, 39421 37019, 37020, 37021, 37022, 37023, 37024, 37025, 37043, 37071, 37089, 37118, 37205, 37206, 37219, 37220, 37234, 37235, 37243, 37244, 37256, 37265, 37270, 37286, 37313, 37314, 37326, 37338, 37355, 37366, 37371, 37393, 37400, 37401, 37424, 37425, 37439, 37444, 37455, 37460, 37495,
AirSense 11 Elite	39004, 39005, 39100, 39123, 39123, 39123, 39123, 39423 , 39481, 39482, 39461, 39106, 39107
AirSense 10 Respond A	37416
AirSense 10 Respond W	37417
AirSense 11 Respond	39317
AirStart™ 10 APAP	37008, 37009, 37202, 37241, 37365
AirStart 10 CPAP	37002, 37004, 37201, 37240, 37364, 37370
Astral™ 100	27001, 27011, 27012, 27021, 27031, 27032, 27051, 27052, 27061, 27062, 27071, 27072, 27081, 27082, 27091 27003, 27013, 27014, 27023, 27033,
Astral 150	27039, 27053, 27054, 27063, 27064, 27068, 27073, 27074, 27083, 27084, 27088, 27093, 27094
Astral External Battery	27918, 27919
Climate Control kit	36875, 36876, 36905, 36915, 36917

ResMed product	Product codes
LifeChoice Activox™ POC	XYC100, XYC100A, XYC100RA, XYC100RB, XYC102R, XYC100B, XYC100B-PRO, XYC100B-SPT, XYC100B-P4L, 200122
Lumis 100 VPAP S	28001, 28002, 28003, 28004, 28005, 28006, 28007, 28008, 28009, 28010, 28011, 28012, 28013, 28302, 28303, 28304, 28322, 28502
Lumis 100 VPAP ST	28101, 28102, 28103, 28104, 28105, 28106, 28107, 28108, 28109, 28306, 38305, 38307, 28505
Lumis 100 VPAP ST-A	28201, 28221, 28313, 28332, 28513
Lumis 150 VPAP ST	28110, 28111, 28112, 28113, 28114, 28115, 28116, 28117, 28118, 28119, 28120, 28121, 28122, 28123, 28124, 28125, 28308, 28309, 28310, 28311, 28312, 28506, 28507, 28508, 28509
Lumis 150 VPAP ST-A	28202, 28203, 28204, 28205, 28206, 28207, 28208, 28216, 28220, 28222, 28223, 28226, 28315, 28316, 28317, 28318, 28320, 28333, 28312, 28513, 28514
Lumis HFT	28401
Mobi	35001, 35003, 35004, 35005, 35006, 35009, 35013
POC External Battery	35052
ResMed Power Station II (RPS II) and ResMed Power Station II-V2 (RPSII-V2)	24921, 24923, 24925, 24926, 24946, 24947, 24959, 24960, 24961, 24962, 24963, 24964, 24981, 24990, 36821, 36822, 36823, 36824, 37342, 37343, 24701, 24702, 24703, 24704, 24705, 24706, 24701, 24704, 24706
S8™ Auto 25	26104, 26105, 26125, 26126
S8 AutoScore II	33044
S8 AutoSet	33105, 33138
S8 AutoSet II	33129, 33141, 33150
S8 AutoSet Spirit	19423, 19426, 33113, 33121, 33122, 33123, 33124, 33125, 33126
S8 AutoSet Spirit II	33130, 33131, 33132, 33133, 33134, 33135, 33136, 33143, 33144, 33145, 33146, 33147, 33148
S8 AutoSet Vantage	19420, 33112
S8 Elite	19413, 19414, 33021, 33026, 33027, 33028, 33031
S8 Compact	19404, 33030
S8 Elite II	33039, 33040, 33041, 33042, 33043, 33045, 33062, 33072, 33073, 33074, 33075 19405, 19408, 19409, 19410, 33001, 33002, 33003, 33007, 33009,
S8 Escape™	33013, 33024, 33036, 33060
S8 Escape II	33048, 33051, 33052, 33053, 33054, 33061, 33064
S8 Lightweight	33032
S8 Lightweight II	33055
S8 Respond	33128, 33137
S9™ AutoSet	36005, 36015, 36025, 36084, 36085, 36086, 36087, 36088, 36089, 36095, 36105, 36115, 36125, 36145, 36180, 36185, 36205, 36225, 36245, 36265, 36285, 36305, 36325, 36345, 36365, 36375, 36385

ResMed product	Product codes
S9 AutoSet 25	36106, 36116, 36126, 36146, 36206, 36226, 36266, 36296, 36366 36100,
S9 AutoSet CS	36110, 36120, 36140, 36200, 36220, 36240, 36260, 36360 36160, 36170
S9 AutoSet CS PaceWave	36211, 36351, 36390 36015, 36065, 36075 36003, 36013, 36023, 36083,
S9 AutoSet CS-A	36093, 36103, 36113, 36123, 36143, 36183,
S9 AutoSet for Her	36203, 36223, 36243, 36263, 36303, 36343, 36363, 36373
S9 Elite	36001, 36011, 36021, 36081, 36091, 36121, 36141, 36201, 36221, 36261, 36301, 36361
S9 Escape	36002, 36012, 36022, 36092, 36122, 36302, 36362 36007, 36017, 36027, 36037, 36047, 36057, 36097, 36367, 36377
S9 Escape Auto	26008, 26009, 26011, 26013, 26928
S9 VPAP™ Adapt	36331
S9 VPAP Adapt SV	36006, 36016, 36026, 36096
S9 VPAP Adapt SV-A	36004, 36014, 36024, 36094, 36114, 36124, 36144, 36204, 36224, 36244, 36264, 36284, 36304 36008, 36018, 36028, 36098, 36108, 36118, 36128,
S9 VPAP Auto	
S9 VPAP S	36148, 36208, 36228, 36248, 36368, 36378
S9 VPAP ST	36009, 36039, 36049, 36059, 36099, 36119, 36129, 36159, 36169, 36209, 36219, 36239, 36339, 36359, 36379
S9 VPAP ST-A	37418 39318
Sleepmate 10 Auto	37419
Sleepmate 11 Auto	24151, 24152, 24153, 24154, 24155, 24156, 24157, 24158, 24159, 24163, 24164, 24166, 24168, 24170 24172 24140, 24141, 24142, 24143, 24144,
Sleepmate 10 Elite	
Stellar™ 100	24145, 24146, 24147, 24149, 24160, 24161, 24162, 24165, 24167, 24169
Stellar 130	26102, 26103, 26106, 26107, 26108, 26109, 26123, 26127
Stellar 150	26112, 26113, 26114, 26115, 26116, 26117, 26118, 26124 26101, 26111, 26121
VPAP IV	26119, 26120
VPAP IV ST	26110, 26122
VPAP Auto 25	
VPAP S	
VPAP ST	



ResMed

Not all products are available in each country/region. Please check with your regional representative for specific country availability. Please contact ResMed Customer Service with any questions.

Australia: +61 (2) 8884 1000, 1800 658 189 (Toll free)
New Zealand: +0800 737 633 (Toll free)
United Kingdom: +44 (1235) 862 997
EU, France: +33 (0) 1 60 183 600
United States: +1 (800) 424 0737 (Toll free)

Kind regards,

Sheila Bruschi
Director, Regulatory Affairs
Sheila.bruschi@resmed.com

		
Stegemannstraße 5 – 7 56068 Koblenz	Telefon: 0261 13388-0 Telefax: 0261 13388-22	
 www.burbach-goetz.de	 info@burbachgoetz.de	